

NYSDOH Environmental Laboratory Approval Program – TO-15 Checklist

LAB ID:	LAB NAME:
DATE:	ASSESSOR NAME:

TO-15 Checklist Determination of VOCs in Air by GC-MS
Method Number:
SOP Number:
Revision Number:
SOP Date:
Personnel records observed:
Data records observed:

NYSDOH Environmental Laboratory Approval Program – TO-15 Checklist

Checklist Category	Ref.	Y	N	NA	ELAP Code	Comment
Required Apparatus and Reagents - Sampling/Concentrator System	7.2.1					TO1000
Does the lab have Electronic Mass Flow Controllers to maintain constant flows of purge, carrier, and sample gas, and to provide an analog output to monitor flow anomalies?	7.2.1.1				TO1001	
Does the autosampler system have a Vacuum Pump , capable of reducing the downstream pressure of the flow controller to provide the pressure differential necessary to maintain controlled flow rates?	7.2.1.2				TO1002	
Are the system's Tubing and Fittings made of Stainless Steel coated with fused silica to minimize active adsorption sites?	7.2.1.3				TO1003	
Does the lab have Stainless Steel Cylinder Pressure Regulators (standard, two-stage cylinder regulators with pressure gauges) for the Helium carrier gas?	7.2.1.4				TO1004	
Does the system include Gas Purifiers to remove organic impurities and moisture from gas streams?	7.2.1.5				TO1005	
Does the system include a Six-port Gas Chromatographic Valve for routing sample and carrier gas flows?	7.2.1.6				TO1006	
Does the system include an appropriate adsorbent trap (Multisorbent Concentrator) packed with an adsorbent having various retentive properties for trace gases?	7.2.1.7				TO1007	
Does the system include a Cryogenic Concentrator which facilitates the separation of compounds on the GC column?	7.2.1.8				TO1008	
Gas Chromatograph / Mass Spectrometric (GC/MS) System Requirements	7.2.2					TO1009
Does the analytical system include a GC that is interfaced to a concentrator and have all required accessories including analytical columns and gases? All GC carrier gas lines must be constructed from stainless steel or copper tubing . Thread sealants made of Non-polytetrafluoroethylene (PTFE) or flow controllers with Buna-N rubber components must not be used .	7.2.2.1				TO1010	
Is the GC equipped with Chromatographic Columns that are recommended for the separation of target nonpolar compounds?	7.2.2.2				TO1011	
Does the system include either a linear quadrupole or ion trap Mass Spectrometer ? It must be capable of scanning from 35 to 300 amu and producing a mass spectrum which meets all the instrument performance acceptance criteria when 50 nag or less of p-bromofluorobenzene (BFB) is analyzed.	7.2.2.3				TO1012	
Is the GC/MS Interface constructed of all-glass, glass-lined, or fused silica-lined materials? Glass and fused silica should be deactivated.	7.2.2.4				TO1013	
Does the Data System (Computer) have software that allows searching any GC/MS data file for ions of a specified mass and plotting such ion abundances versus time or scan number? This type of plot is defined as a Selected Ion	7.2.2.5				TO1014	

NYSDOH Environmental Laboratory Approval Program – TO-15 Checklist

Checklist Category	Ref.	Y	N	NA	ELAP Code	Comment
Current Profile (SICP). Software must also be available that allows integrating the abundance in any SICP between specified time or scan number limits. Also, software must be available that allows for the comparison of sample spectra with reference library spectra.						
Is the Off-line Data Storage Device capable of rapid recording and retrieval of data, and is it also suitable for long-term, off-line data storage?	7.2.2.6				TO1015	
Calibration System and Manifold Apparatus	7.3	TO1016				
Does the system have a Calibration Manifold used for the preparation of gaseous standards (either single compounds or multicomponent gases)?	7.3.1				TO1017	
Does the system have a Humidifier containing HPLC grade deionized water? Ideally, the HPLC water should be boiled prior to use in order to ensure VOC-free water.	7.3.2				TO1018	
Does the system have Electronic Mass Flow Controllers For the calibration gases?	7.3.3				TO1019	
Does the system have Filters inside the manifold for the collection of particulates? (Usually these are 47-mm, made of Teflon®)	7.3.4				TO1020	
Reagents	7.4	TO1021				
Does the lab have Neat Gases or Manufacturer-Certified Gas Mixtures ?	7.4.1				TO1022	
Does the lab use Ultra-high purity grade Helium for the carrier gas in the GC?	7.4.2				TO1023	
Does the lab have Liquid Nitrogen or Liquid Carbon Dioxide used for cooling the secondary trap?	7.4.3				TO1024	
Does the lab use HPLC grade Deionized Water for the humidifier?	7.4.4				TO1025	
Canister Cleaning and Certification	8.4.1	TO1026				
Are all canisters clean and free of any contaminants before sample collection? Review the cleaning & certification documentation related to each canister by ID tag.	8.4.1.1				TO1027	
Are all canisters leak tested by pressurizing them to approximately 206 kPa (+30 psig) with humid zero air or zero N ₂ ? Alternatively, are the canisters evacuated (to -30psig vacuum)? The initial pressure (or vacuum) is measured, and after 24 hours, the final pressure (or vacuum) is checked. To be acceptable, the reading should not vary more than ± 13.8 kPa (± 2 psig) over the 24 hour period.	8.4.1.2				TO1028	
Has any canister that has not tested clean at less than 0.2 ppbv of targeted VOCs not been used, or have the contaminants been noted in the cleaning/certification logbook?	8.4.1.6				TO1029	
As an option to the humid zero air cleaning procedures, are the canisters heated in an isothermal oven between 100°-105°C during evacuation of the canister? At the end of the heating / evacuation cycle, the canisters are pressurized with humid zero air (or gaseous nitrogen) and analyzed by a	8.4.1.8				TO1030	

NYSDOH Environmental Laboratory Approval Program – TO-15 Checklist

Checklist Category	Ref.	Y	N	NA	ELAP Code	Comment
GC/MS system after a minimum of 12 hrs of "aging."						
Cleaning of the Sampling System Components	8.4.2	TO1031				
Are sampling system components disassembled and cleaned before the sampler is assembled?	8.4.2.1				TO1032	
Are the parts then rinsed with HPLC grade deionized water and dried in a vacuum oven at 100°C for 12 to 24 hours?	8.4.2.2				TO1033	
Once the sampler is assembled, is the entire system purged with humid zero air for 24 hours?	8.4.2.3				To1034	
Zero Air Certification	8.4.3	TO1035				
Is the cleanliness of the sampling system determined by testing the sampler with humid zero air or humidified zero nitrogen? (N ₂ is 100% free of VOCs and need not be passed through a chemical trap).	8.4.3.1				TO1036	
If the sampler passes the humid zero air test, is it then tested with humid calibration gas standards containing selected VOCs at concentration levels expected in field sampling (e.g., 0.5 to 2 ppbv)?	8.4.3.3				To1037	
Not all the target gases on the Title III list are available or compatible with compressed gas standards. In these cases, is sampler certification approached by a different means?	8.4.5.1				TO1038	
Preparation of Standards	9.2	TO1039				
Are the standard mixtures of target gases in high pressure cylinders certified traceable to a NIST or EPA Certified Reference Materials?	9.2.1.1				TO1040	
Are the "neat" standards that are used for making trace gas standards of high purity? Generally, a purity of 98% or better is commercially available	9.2.1.2				TO1041	
Does the lab retain the traceability certificates for the cylinders containing primary stock standards? (These are usually at concentrations from 100 ppbv to 1000 ppbv)	9.2.1.3				TO1042	
Does the lab have an Instrument Performance Check Standard? This is a prepared standard solution of BFB in Nitrogen or humidified zero air, at a concentration which will allow collection of 50 ng of BFB or less under the optimized concentration parameters.	9.2.2.1				TO1043	
Does the lab prepare five working calibration standards at the 2, 5, 10, 20, and 50 ppbv level for each component? A separate cal curve ranging from approximately 0.5 ppbv to 2 ppbv may be appropriate for labs with low-level capability.	9.2.2.2				TO1044	
Does the lab have an Internal Standard Spiking Mixture: A.) Is an internal spiking mixture containing bromochloromethane, chlorobenzene-d, and 1,4-difluorobenzene added to all samples and standards? B.) Is the Int Std introduced into the trap during the collection time for all calibration, blank, and sample analyses? Alternatively, ISTD containing BFB as the fourth component allows easy evaluation of BFB tune in any sample or blank, plus it allows using BFB as a retention time locking compd. C.) Is the volume of internal std spiking mixture added for	9.2.2.3				TO1045 TO1046	

NYSDOH Environmental Laboratory Approval Program – TO-15 Checklist

Checklist Category	Ref.	Y	N	NA	ELAP Code	Comment
each analysis the same from run to run?					TO1047	
Are the Standards prepared by dynamic dilution of the calibration stock standards with zero Nitrogen or humidified zero air, using mass flow controllers and a calibration manifold?	9.2.3.1				TO1048	
Does the lab have a Standard Preparation Procedure for use in High Pressure Cylinders? Usually, the standards are prepared as gas mixtures in a Dynamic Diluter. Gases are mixed in the diluter and then transferred to a cylinder to approx 20 psia or 5 psig. Alternatively, predetermined amounts of each neat standard compound are measured using a microliter or gastight syringe and injected into the cylinder.	9.2.5				TO1050	
Does the lab have a procedure for the Storage of Working Standards? Working standards prepared in canisters may be stored for 30 days.	9.2.8				TO1052	
Gas Chromatograph / Mass Spectrometric (GC/MS) Operating Conditions	10	TO1053				
Does the lab use the following Analytical Sequence? <ul style="list-style-type: none"> • Perform instrument performance check using BFB. • Initiate multi-point calibration or daily calibration checks. • Perform a laboratory method blank. • Complete this sequence for analysis of 20 field samples or each 24 hr period. 	10.3				TO1054	
Does the lab use a BFB Instrument Performance Check (IPC)? The instrument performance check solution must be analyzed initially and once per 24-hour time period of operation	10.4				TO1055	
Does the MS Tune pass the abundance criteria for BFB , and retune the Mass Spec if BFB acceptance criteria aren't met?	10.4.4				TO1056	
Are the Results of the BFB tuning recorded and maintained as part of the instrumentation log?	10.4.6				TO1057	
Prior to the analysis of samples & blanks, but after the IPC criteria have been met, is each GC/MS system calibrated at 5 concentrations that span the monitoring range of interest (or is a daily CCV checked)? One of the calibration points from the initial calibration curve must be at the same concentration as the daily calibration standard (e.g., a midrange 10 ppbv).	10.5				TO1058	
Is each GC/MS system recalibrated following corrective action (e.g., ion source cleaning or repair, column replacement, etc.) which may change or affect the initial cal criteria, or if the daily calibration acceptance criteria have not been met?	10.5.2				TO1059	
Concentration of cal standards: A minimum of five concentration levels are needed to determine the instrument sensitivity and linearity. Is one of the calibration levels near the detection level for the compounds of interest?	10.5.3				TO1060	
Are Initial calibration acceptance criteria met (+/-30% D)	10.5.6.2				TO1061	

NYSDOH Environmental Laboratory Approval Program – TO-15 Checklist

Checklist Category	Ref.	Y	N	NA	ELAP Code	Comment
before any field samples, blanks, or performance evaluation (PE) samples are analyzed?						
Prior to analyzing samples and blanks but after tuning criteria are met, is the ICAL checked by analyzing a daily cal std to ensure that the instrument is under control? The daily calibration standard should contain all the target compounds.	10.6				TO1062	
Is a check of the calibration curve performed once every 24 hours of run time on the GC/MS system?	10.6.2				TO1063	
Has the mid-level cal standard (10 ppbv) analyzed by the GC/MS system met the tuning and mass calibration Criteria?	10.6.3				TO1064	
Is the Relative Response Factor (RRF) calculated for each target Compound?	10.6.4.1				TO1065	
Percent Difference (%D): The percent difference in the daily RRF (24-hour) must be compared to the mean RRF in the most recent initial calibration. Does the lab calculate %D for each target compound?	10.6.4.2				TO1066	
Technical Acceptance Criteria: The %D for each target compound in a daily calibration sequence must be within ± 30% to proceed with analysis. Does the lab maintain a control chart showing %D values?	10.6.5				TO1067	
Blank Analyses: To monitor for possible laboratory contamination, does the laboratory analyze a Method Blank at least once in a 24-hour analytical sequence? All steps in the analytical procedure are performed on the blank using all items that would be used for a sample analysis.	10.7				TO1068	
Does the laboratory analyze a method blank after the calibration standard(s) and before any samples are analyzed?	10.7.2				TO1069	
Technical Acceptance Criteria for daily runs: <ul style="list-style-type: none"> • Is a blank canister analyzed daily? • Is the area response for each internal standard in the blank within ± 40% of the mean area response of the Internal Std in the most recent valid calibration? Also, is the retention time for each of the internal standards within ±0.33 minutes between the blank and the most recent valid calibration? • Does the blank NOT contain any target analyte at a concentration greater than its quantitation level (three times the MDL)? Also, does the blank NOT contain additional compounds with elution characteristics and mass spectral features that would interfere with identification and measurement of a method analyte? 	10.7.5				TO1070 TO1071 TO1072	
Corrective Action: If the blanks do not meet the technical acceptance criteria, the analyst should consider the analytical system to be out of control. If a blank is contaminated with an analyte, and that analyte is also found in associated samples, are those sample results "flagged" as possibly contaminated?	10.7.6				TO1073	
If time remains in the 24-hour period in which an initial cal is	10.8.2				TO1075	

NYSDOH Environmental Laboratory Approval Program – TO-15 Checklist

Checklist Category	Ref.	Y	N	NA	ELAP Code	Comment
performed, samples may be analyzed without analysis of a daily cal standard. If time does not remain in the 24-hr period since the injection of the instrument performance check std during which the initial calibration was performed, are both the IPC standard and the daily calibration standard analyzed before sample analysis begins?						
Secondary ion quantitation is allowed only when there are sample matrix interferences with the primary ion. If secondary ion quantitation is performed, is the reason documented in the laboratory record book?	10.8.3.9				TO1076	
Are Field samples analyzed only after the GC/MS system meets the BFB tuning, initial calibration, and continuing calibration technical acceptance criteria?	10.8.5.1				TO1077	
Are the field samples analyzed along with a lab Method Blank that met the blank technical acceptance criteria?	10.8.5.2				TO1078	
Are all of the target analyte peaks in the samples within the initial calibration range?	10.8.5.3				TO1079	
Is the retention time for each internal standard within ±0.33 minutes of the retention time of the internal standard in the most recent valid calibration?	10.8.5.4				TO1080	
If the on-column concentration of any compound in any sample exceeds the initial calibration range , is an aliquot of the original sample diluted and reanalyzed?	10.8.6				TO1081	
Internal standard responses and retention times must be evaluated after data acquisition. If the retention time for any internal standard changes by more than 20 sec from the latest daily (24-hour) calibration standard (or mean retention time over the initial calibration range), is the GC/MS system inspected and corrections made?	10.8.6.1				TO1082	
If the area response for any internal standard changes by more than ± 40% between the sample and the most recent valid calibration, the GC/MS system must be inspected for malfunction and corrections made as appropriate. When corrections are made, are those samples reanalyzed?	10.8.6.2				TO1083	
If, after reanalysis, the area responses or the RTs for all internal standards are inside the control limits, then the problem with the first analysis is considered to have been within the control of the Laboratory. Does the lab submit only data from the analysis with SICPs within the limits? This is considered the initial analysis and should be reported as such on all data deliverables.	10.8.6.3				TO1084	
Requirements for Demonstrating Method Acceptability for VOC Analysis from Canisters	11.	TO1085				
Does the lab meet the three performance criteria for an analytical system to qualify under Compendium Method TO-15? These criteria are: <ul style="list-style-type: none"> • an MDL of no greater than 0.5 ppbv, • replicate precision within 25%, and • audit accuracy within 30% for concentrations expected 	11.1.1				TO1086	

NYSDOH Environmental Laboratory Approval Program – TO-15 Checklist

Checklist Category	Ref.	Y	N	NA	ELAP Code	Comment
in contaminated ambient air (e.g. 0.5 to 25 ppbv).						
Is the procedure chosen to define the method detection limit the one specified in the Code of Federal Regulations (40 CFR 136 Appendix B)?	11.2				TO1087	
Is the measurement of replicate precision calculated according to the following procedure: The absolute value of the difference between replicate measurements of the sample divided by the average value, and expressed as a percentage (%RPD)?	11.3				TO1088	

APPENDIX B. COMMENT ON CANISTER CLEANING PROCEDURES

The canister cleaning procedures given in Section 8.4 require that canister pressure be reduced to <0.05mm Hg before the cleaning process is complete. Depending on the vacuum system design (diameter of connecting tubing, valve restrictions, etc.) and the placement of the vacuum gauge, the achievement of this value may take several hours. In any case, the pressure gauge should be placed near the canisters to determine pressure. The objective of requiring a low pressure evacuation during canister cleaning is to reduce contaminants. If canisters can be routinely certified (<0.2 ppbv for target compounds) while using a higher vacuum, then this criteria can be relaxed. However, the ultimate vacuum achieved during cleaning should always be <0.2mm Hg. Canister cleaning as described in Section 8.4 requires components with special features. The vacuum gauge must be capable of measuring 0.05mm Hg with less than a 20% error. The vacuum pump used for evacuating the canister must be noncontaminating while being capable of achieving the 0.05 mm Hg vacuum as monitored near the canisters. Thermoelectric vacuum gauges and turbomolecular drag pumps are typically being used for these two components.

An alternate to achieving the canister certification requirement of <0.2 ppbv for all target compounds is the criteria used in Compendium Method TO-12 that the total carbon count be <10ppbC. This check is less expensive and typically more exacting than the current certification requirement and can be used if proven to be equivalent to the original requirement. This equivalency must be established by comparing the total non-methane organic carbon (TNMOC) expressed in ppbC to the requirement that individual target compounds be <0.2 ppbv for a series of analytical runs.