

## Blood pH and Gases

<b>Blood pH and Gases</b>	
<b>Standard</b>	<b>Guidance</b>
<p><b>Blood pH and Gases Standard of Practice 1 (BP S1): Quality Control Requirements</b></p> <p>Unless an Individualized Quality Control Plan (IQCP) is established according to <a href="#">Quality Control Standards of Practice 2, 3 and 4</a>, for all blood gas instruments the laboratory must:</p> <ul style="list-style-type: none"><li>a) test three (3) levels of controls on each day of testing, and a minimum of one (1) control each eight (8) hours of testing; and</li><li>b) test a control or calibrator with each patient unless the blood gas instrument internally verifies the calibration every thirty (30) minutes at a minimum.</li></ul>	